

Mosaic 9 Integrations

Mosaic software helps to digitalise your lab by offering drivers and support for tight 'driven' integrations with an extensive range of instruments and devices from multiple vendors. Auto-validated file integration is possible with any lab automation

 **AltemisLab**

analytikjena
An Endress+Hauser Company

 **AZENTA**
LIFE SCIENCES

 **BECKMAN
COULTER** *Life Sciences*

BioMicroLab

 **BIOSERO**

 **BIOVIA**

 **CDD.VAULT®**
Complexity Simplified

 **ChemAxon**

dotmatics
knowledge solutions

Genedata
solutions in silico 

HAMILTON

 **HighRes®**
biosolutions

METTLER TOLEDO

 **MICRONIC**
INNOVATION IN SAMPLE STORAGE

revvity
signals

SARTORIUS

 **sptlabtech**

 **TECAN.**

ThermoFisher
SCIENTIFIC

Ziath

MOSAIC SOFTWARE SUPPORTS THE FOLLOWING:

Automated Stores

- Azenta Life Sciences SampleStore™
- Azenta Life Sciences BioStore™
- Azenta Life Sciences LN2-Based BioStore™ -80 and -190
- Azenta Life Sciences LN2-Based Cryo Store Pico™
- Azenta Life Sciences Matrical MiniStore
- Hamilton Storage BIOS
- Hamilton Storage SAM HD
- Hamilton Storage Verso®
- Hamilton Storage Verso® Q20
- SPT Labtech compACT
- SPT Labtech comPOUND®
- SPT Labtech arktic®

Manual Stores

- Kardex Shuttle

Biological Registration

- Dassault Systemes BIOVIA™ Registration
- Dotmatics BioRegister
- Genedata Biologics
- Watson LIMS

Screening Data Analysis

- Dotmatics Ultra Screening
- Genedata Screener
- IDBS ActivityBase
- Rewity Signals VitroVivo

Small Molecule Registration

- CDD Vault Registration System
- ChemAxon Chemical Registration
- Dassault Systemes BIOVIA™ Chemical Registration
- Deltasoft ChemCart
- Dotmatics Register
- IDBS ActivityBase
- Scilligence RegMol

Other Software

- Certara D360
- Dassault Systemes BIOVIA™ Workbook
- Dotmatics Browser

Dry Powder Handling

- Mettler-Toledo balances
- Mettler-Toledo FlexiWeigh
- Sartorius balances

2D Tube Sorting

- SPT Labtech BioMicroLab XL9
- SPT Labtech BioMicroLab XL20
- SPT Labtech BioMicroLab XL100
- SPT Labtech BioMicroLab XL200
- Ziath Mohawk



Liquid Handling

- Agilent BioCel
- Agilent Bravo
- Beckman Coulter Life Sciences Biomek FX
- Beckman Coulter Life Sciences Biomek NX
- Beckman Coulter Life Sciences Biomek i-Series
- Beckman Coulter Life Sciences Echo 525 & 650 (Standalone)
- Beckman Coulter Life Sciences Access Laboratory Workstation
- Beckman Coulter Life Sciences Access Dual Robot System
- Biosero Green Button Go® Scheduler
- Hamilton STAR
- Hamilton STARlet
- Hamilton STARplus
- Hamilton Vantage
- HighRes Biosolutions ACell
- HighRes Biosolutions CoLab Flex
- HighRes Biosolutions MicroStar
- HighRes Biosolutions NanoCell
- Matrical AutoMap
- Perkin Elmer MultiPROBE II
- Tecan Freedom EVO®
- Tecan Fluent®
- SPT Labtech mosquito®
- Zinsser Lissy

Other

- Any “keyboard wedge” 2D barcode reader
- Any label printer via BarTender labelling
- Azenta Life Sciences Tube Auditor
- CyBio CyBi-Print labeller
- LogistiVIEW Digital Workforce Platform
- SPT Labtech BioMicroLab VC100
- SPT Labtech BioMicroLab VC384
- Vuzix M300, M300XL & M400 smart glasses

2D Rack Scanning

- AltemisLab AlteRead Rack Scanner
- Azenta Life Sciences Scanner-Based Barcode Reader
- Azenta Life Sciences Camera-Based Reader for SBS Racks
- Azenta Life Sciences Camera-Based Reader for Acoustic Tubes
- Biosero BioRead A6 High Speed Single Rack Scanner
- Biosero BioRead AV6 Deep Focus Single Rack Scanner
- Biosero BioRead Express Rack Reader
- Biosero BioRead Multirack Scanner
- Biosero DataPaq Cube Scanner
- Hamilton Storage ID Reader
- Micronic Tracxer Code Reader
- Micronic DR700 Code Reader
- Micronic DR710 Code Reader
- Micronic DR505 Code Reader
- SPT Labtech BioMicroLab Scan + Barcode Scanner
- SPT Labtech BioMicroLab Scan M Barcode Scanner
- SPT Labtech BioMicroLab Scan HS Barcode Scanner
- SPT Labtech BioMicroLab Scan ES Barcode Scanner
- SPT Labtech BioMicroLab Scan C Barcode Scanner
- Thermo Scientific VisionMate™ HS 2D Barcode Reader
- Thermo Scientific VisionMate™ HSX 2D Barcode Reader
- Zebra TC series touch computer
- Ziath DataPaq™ Cube Rack Reader
- Ziath DataPaq™ Deep Focus Rack Reader
- Ziath DataPaq™ Express Rack Reader
- Ziath DataPaq™ High Speed Scanner
- Ziath DataPaq™ Multirack Scanner
- Ziath DataPaq™ Single Tube Scanner





Different device?

Almost any lab equipment can be integrated with Mosaic software using our validated file exchange or JFM Application Programming Interfaces.

If you don't see the device you want in this document, please talk to a Titian representative.

Thank you for your interest in this document



White papers



Application notes



Product brochures



Video content



Webinar recordings

You can find a collection of further related material in the resources section on our website.

info@titian.co.uk

www.titian.co.uk

UK Tel:
+44 20 7367 6869

USA Tel:
+1 508 366 2234